

P O W E R C H A I R



Seat Colors
Gray Vinyl

Body Colors
Candy Apple Red
Blue

Roll



New and Improved Controller!

golden COMPASSTM SPORT

The 2015 Golden Compass Sport features the new Dynamic Linx controller! The new Linx offers a better ride experience using patented Dynamic Load Compensation technology. This new technology improves the terrain-climbing ability of the Compass Sport's patented articulating chassis. Learning to drive the Compass Sport is easier and a more natural experience than before!



Rise.

Roll.

Ride.



YOUR LIFE IN MOTION

Golden Compass™ Sport Specifications

Model	GP605CC	GP605SS
Model Name	Compass Sport	Compass Sport
Medicare Code	K0823	K0822
Maximum Incline	6°	6°
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Center	Center
Maximum Speed	4.5 mph	4.5 mph
Operating Range ¹	19.5 miles ³	19.5 miles ³
Ground Clearance	3.25"	3.25"
Turning Radius	19.5"	19.5"
Type Batteries	2-U1	2-U1
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	Dynamic Linx 40 Amp	Dynamic Linx 40 Amp
Dimensions:		
Length	39"	39"
Width (outside of tires)	24"	24"
Height (ground to top of back of seat)	38" - 40"	40" - 42"
Height (ground to top of back of headrest)	45" - 50.5"	47" - 52.5"
Ground to Top of Seat	21" - 23"	19" - 21"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Weight of Unit: (assembled)		
Base	92 lbs	92 lbs
Seat (18 x 18 Solid Pan Seat)		40 lbs
Seat (18 x 18 Captain's Chair Seat)	46 lbs	
Batteries	U1 - 22 lbs	U1 - 22 lbs
Tires:		
Flat Free	Flat Free	Flat Free
Center Drive Wheels (2)	10"	10"
Casters (4)	6"	6"
Standard Seat:		
Captain's Seat	Captain's Seat	Captain's Van/Pan
Back Height - no headrest	19"	23"
Back Height with headrest	26.5" - 29.5"	30.5" - 34.5"
Width x Depth (Adjustable Depth)	18" x 18"	18" x 18"
Color	Gray Vinyl	Coffee/Sand, Charcoal/Gray
Footrest:		
Yes	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
Length adjustable (Option)	Yes	Yes
Warranty:		
Frame	Lifetime	Lifetime
Drive Train	13 Months	13 Months
Electronics ²	13 Months	13 Months

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

³ Range listed tested at 300 lbs.

Armrests are height, width and angle adjustable



2 gorgeous colors



Candy Apple Red



Blue

Available Accessories:

- Crutch Holder Tube
- Cane Holder Tube
- Oxygen Tank Holder
- Quad Cane Holder
- Walker Holder
- Swing-away Desk Arm
- Swing-away Footrests
- Stump Support
- Elevating Leg Rests
- Power Chair Cover
- Rear Basket
- Pack N' Go
- Van/pan Seat
- EZ Lock System